

Aqua Creek Products

The Total Solution for Pool and Spa Access



ADA Solution For Existing Anchors

Independent third party verified to be ADA Compliant, The Scout™ is specifically designed to be used in retro-fit applications. By using one of our multiple adapter sleeves, The Scout can fit a wide variety of pre-existing lift anchors.

The Scout incorporates a sturdy stainless steel structure.




Item# F-800SC-DER-D

U.S. Patent No.
D507,769 S

OPTIONS

- **Upgrade Pack - Item #F-800SCAP**
Includes battery, headrest, chest strap & cover
- **Headrest - Item #F-722HRR**
- **Chest strap - Item #F-423CS**
- **Pull-out leg rest - Item #F-105LAR**
- **Lockable cover - Item #F-450SC**
- **Anchor, standard (additional) - Item #F-808SA**
- **Transport Cart - Item #F-814SCTC**

FEATURES

- Independently Verified to ADA Standards
- 350 lb* (158.5 kg) weight capacity
- Set back range 12" to 36" (30.5 to 91 cm)
- Adaptable to most existing 6" deep anchor sockets
- 360° powered rotation in either direction
- Rechargeable battery with wall mount charger
- Limited lifetime structural warranty
- 5 year pro-rated electronics warranty
- CE version available 

**300 lb (136 kg) capacity when used with 1.9" anchors with plastic sleeves*



PROUDLY SOLD BY: