

**NEW**

# Transforming sleep.

## Forma™

### FULL FACE MASK

The all-new Forma™ Full Face Mask forms a superior seal and provides even more comfort with a greater range of movement through the active contouring of its new FlexiFoam™ Cushion. A new T-Piece adds stability to the stylized Mask Base, and works in harmony with the highly contoured Silicone Seal and Under-Chin Design to enhance sleep performance. With its ease of use and superior conforming seal, the Forma™ sets a whole new standard.


Transform your patient's sleep experience today, with the Fisher & Paykel Healthcare Forma™ Full Face Mask.

NEW:: FlexiFoam™  
& T-Piece Mask Base




**Fisher & Paykel**  
HEALTHCARE

# Full Face Delivery – *Forma*<sup>™</sup>




The *Forma*<sup>™</sup> Full Face Mask offers an improved seal and a more comfortable fit due to the interaction between the new FlexiFoam<sup>™</sup> Cushion and the Under-Chin Design. This alternative to nasal or oral delivery comes fully assembled and ready for use.

## UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE




### FlexiFit<sup>™</sup> Technology Auto-contouring fit

The FlexiFoam<sup>™</sup> Cushion using FlexiFit<sup>™</sup> Technology automatically adjusts to a wide range of nasal and facial contours with simple one-step ease of fitting and no complicated adjustments.




### New:: FlexiFoam<sup>™</sup> Cushion Active-contouring cushion

The new FlexiFoam<sup>™</sup> Cushion provides greater comfort with the unique design of the soft, light and supportive foam. The FlexiFoam<sup>™</sup> features the advanced active contouring that conforms naturally to your patient's face.




### Glider<sup>™</sup> Strap with Easy-Release Cord

The sliding Glider<sup>™</sup> Strap with Easy-Release Cord provides freedom of movement and one-step removal. The mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks. Simply pull the Easy-Release Cord for one-step removal.




### Under-Chin Design with Crown Strap

An Under-Chin fit provides extra stability and maintains the seal, even with movement of the jaw during the night. The Crown Strap works in conjunction with the Under-Chin Design to eliminate the need for a chin strap.



### New:: T-Piece Mask Base Optimal seal

With no complicated adjustments, the T-Piece provides added stability and complements the Under-Chin Design to allow FlexiFit<sup>™</sup> to auto-adjust for hassle-free fitting.



### New:: Advanced Air Diffuser Quiet usage

*Maintenance Free.* This efficient Air Diffuser system prevents air from blowing onto the sleeping partner and decreases noise that could potentially disrupt sleep.



### Frosted Silicone Seal Optimizing seal

The seal technology has a soft, smooth feel to slide onto the face and its frosted texture optimizes the mask seal.

## CONTENTS

- Complete mask with headgear assembled and ready for use
- Four sizes: Small, Medium, Large and X-Large, sold separately
- User-friendly Fitting Guide

REORDER CODE: 400470 Small  
400471 Medium  
400472 Large  
400473 X-Large

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER & PAYKEL HEALTHCARE REPRESENTATIVE:

**Manufacturer** Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348, Panmure, Auckland 1741, New Zealand  
Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web: www.fphcare.com

**Australia** Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **China** Tel: +86 20 3205 3486 Fax: +86 20 3205 2132  
**France** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99 **India** Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 **Irish Republic** Tel: 1800 409 011  
**Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **Middle East** Tel: +9 71 4 360 0849 Fax: +9 71 4 368 8110 **Northern Ireland** Tel: 0800 132 189  
**Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK (EU Authorised Representative)**  
Fisher & Paykel Healthcare UK Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: +44 1628 626 136 Fax: +44 1628 626 146  
**USA** Tel: 1 800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001